

## UF SERIES HOLLOW FIBER POLYSULFONE MEDIA FILTER CARTRIDGE

- Composed of validatable hollow fiber polysulfone media with polypropylene hardware
- These absolute rated filters offer excellent bacteria, particle and pyrogen removal for critical pharmaceutical and life science applications, at economical costs
- Code 0, Code 5, and Code 7 configurations are available in 10, 20, 30 and 40-in lengths
- Meet current USP Class VI biological test for plastics, with a validation certificate available
- Contain no surfactants or coatings to interfere with your application, and thermo-plastic construction process minimizes extractables
- Hollow fiber construction ensures strength, while exhibiting low pressure drop and exceptionally high flow rates



## **CARTRIDGE SPECIFICATIONS**

| Dimensions                             |  |
|--|--|
| Diameter (OD)                          | 2.9" (7.4 cm)  |
| Length (in)                            | 10, 20, 30, 40   |
| Materials                              |  |
| Cartridge                              | Polypropylene caps and cage, polyurethane end sealant  |
| Media                                  | Polysulfone fibers                                     |
| Seals                                  | Silicone o-rings                                       |
| Operating Parameters                   |  |
| Maximum Temperature                    | 185°F (85°C) @ 10 psid (.7 bar)                        |
| Maximum Differential Pressure          | 30 psi (2 bar) @ 25°C                                  |
| Filtration Rating                      | 100% Absolute  |
| Toxicity                               | Non-toxic by USP Class VI Biological Test for Plastics |
| Validation                             | Diffusive flow test methods                            |
| Recommened Water Flow Rates            |  |
| 0.05 micron                            | 2.0 gpm/10" length (7.6 lpm/10" length)                |
| 0.1 and 0.2 micron                     | 3.0 gpm/10" length (11.4 lpm/10" length)               |
| Sanitizing Agents                      |  |
| Chlorine, hydrogen peroxide, autoclave | d for 20 minutes at 121°C under no load conditions     |

## **ORDERING INFORMATION**

| Catalog Number and Description |   |  |
|--------------------------------|---|--|
| FCUF                           | Filter Cartridge Ultra Filter   |  |
| X                              | Cartridge Code: <b>0</b> = 2-222 o-rings <b>5</b> = 2-222 o-rings w/spear end <b>1</b> = double open end w/ internal o-ring <b>7</b> = 2-226 o-rings w/locking tabs/spear end |  |
| XX                             | Length (in): <b>10</b> = 10, <b>20</b> = 20, <b>30</b> = 30, <b>40</b> = 40   |  |
| XX(X)                          | Micron Rating: <b>S05</b> = 0.05 $\mu$ m, <b>S1</b> = 0.1 $\mu$ m, <b>S2</b> = 0.2 $\mu$ m  |  |

To figure your order number, replace the X with one of the numbered or lettered options beside it. Note: Not all part number combinations are available; consult Technical Support for assistance.

24-Hour Customer Service: +1 (800) 466-7873 www.evoqua.com

Water Flow Rate/10" Cartridge (GPM)

FLOW RATE vs. INITIAL CLEAN PRESSURE DROP

All information presented herein is believed reliable and in accordance with accepted engineering practices. Evoqua makes no warranties as to the completeness of this information. Users are responsible for evaluating individual product suitability for specific applications. Evoqua assumes no liability whatsoever for any special, indirect or consequential damages arising from the sale, resale or misuse of its products.